




Uromic QuickStep

Multifunctional urodynamic system provides "Best in Class" specification.

- Compact 8 channel system for complex functional diagnostics.
- Simple hydraulic height adjustment of the keyboard and monitor enables seated or standing position.
- Small footprint minimises use of space in clinic and storage.
- Touchscreen monitor and infection control friendly keyboard.
- Unique wireless uroflowmeter connection at 866/915 MHz ensures the stability and no signal dropout.
- Further options for videourodynamics available.
- Accurate 7 roller pump.



↑
Hydraulic pedal for easy height adjustment and 4 braked wheels for safety.



↑
Option: inbuilt ultrasound module for complete pelvic floor examination.



↑
Adapter for wireless connection at 866/915 MHz.