

High-Performance, All-in-One Ultrasound System



Excellent Performance in a Compact, Portable Design

- High Resolution Imaging
- Advanced Imaging Technologies
- Fully Customizable Touch Screen
- Versatile Applications and Easy-to-Use Features



Quantum
Technology

Maximize your Imaging Performance with Quantum Plus Technology

SonixOne is an easy-to-use, fully featured, compact ultrasound system incorporating Quantum Plus Technology for excellent image quality.

High Resolution Imaging with Advanced Technologies

- Efficient and fast performance
- High resolution capacitive 19" touch screen monitor with advanced zero bezel design
- Sensitive color Doppler with excellent spatial resolution and sensitivity
- Motion Compensated Angular Compound Imaging (MACI) ensures high quality real-time images
- Multiple imaging modes including Spatial Compounding, Trapezoidal, and Panoramic Imaging

Versatile Applications

SonixOne is a flexible system ideal for many applications including:

- Oocyte retrievals for IVF
- Biopsy guidance in interventional suites
- Small clinics with limited budgets
- Fine needle aspiration
- Specialty applications such as podiatry and MSK
- Secondary systems where space is at a premium

Innovative Needle Guidance

- SonixShine needle enhancement technology illuminates needle in-plane with the touch of a button
- Ideal for guiding needle placement during interventional procedures

Portable, Flexible and Optimized Workflow


- Mount on a wall, ceiling or height-adjustable cart using universal VESA support or position on any flat, stable surface
- Fully programmable touch screen for easy customization
- SonixHub streamlines credentialing, billing and quality assurance processes




Save Space Without Compromising Performance
SonixOne can be placed on a stable surface and has a convenient handle for transport



Customizable Touch Screen
Equipped with a configurable touch screen, SonixOne is completely customizable to show only the controls that are relevant to the exam at hand



SonixShine Technology
Illuminate the needle at the touch of a button




SonixHub - Credentialing and Billing Made Easier
Reduce the time, cost and paperwork in your busy practice

Compact with Excellent Performance and Advanced Features



Gallbladder polyp




Abdominal aorta pulsed wave Doppler



Angiomyolipoma at lower pole of right kidney



Thyroid right transverse



Elastography thyroid nodule



Thyroid cystic mass with color Doppler



Ovarian follicles



4 chamber heart - 30 week fetus




Uterus transverse



Calcification left subscapularis



*Transverse right kidney with renal vein.
Gallbladder fold in transverse plane*















Shoulder supraspinatus tendon

Innovative Solutions for Life

Analogic Corporation creates innovative technology to improve the health and ensure the safety of people around the world. We are committed to providing ultrasound solutions under the Ultrasonix brand name that advance medicine and save lives.

High-Performance, Lightweight Transducers

Curved C5-2/60	Curved C7-3/50	Microconvex MC9-4/12	Endocavity Microconvex EC9-5/10	Transrectal Biplane Curved BPC8-4/10	Transrectal Biplane Curved BPL9-5/55
5-2 MHz	7-3 MHz	9-4 MHz	9-5 MHz	8-4 MHz	9-5 MHz
					
Linear L9-4/38	Linear L14-5/38	Linear L14-5W/60	Linear L40-8/12	Linear Hockey Stick HST15-8/20	Phased Array SA4-2/24
9-4 MHz	14-5 MHz	14-5 MHz	40-8 MHz	15-8 MHz	3-2 MHz
					



 Ultrasonix Medical Corporation, 130-4311 Viking Way, Richmond, BC V6V 2K9, Canada. T +1 604-279-8550 F +1 604-279-8559



powered by analogic 

BK Ultrasound
8 Centennial Drive
Peabody
MA 01960
USA
T +1 978-326-4000
bkultrasound.com

Sales and Service USA
BK Ultrasound
8 Centennial Drive
Peabody
MA 01960
USA
T +1 978-326-1300
F +1 978-326-1399
bkultrasound.com

Europe and Rest of World
BK Ultrasound
Mileparken 34
2730 Herlev
Denmark
T +45 4452 8100
F +45 4452 8199
bkultrasound.com

Canada
BK Ultrasound
#130-4311 Viking Way
Richmond BC
Canada V6V2K9
T +1-866-437-9508
F +1-604-279-8559
bkultrasound.com

Asia
Analogic Medical Equipment
(Shanghai) Co., Ltd.
1377, Lan Dian Road,
Pu Dong New District,
Shanghai, China 201132
bkultrasound.com